

Mouse IgG, Protein A Purified, pAb

Catalog No.	CSI20204A CSI20204B	Quantity:	10 mg 50 mg
Description:	Purified from normal serum by immobilized Protein A.		
Concentration:	14.1 mg/mL		
Volume:	3.55 mL		
Purity:	IgG fraction, >95% by SDS-PAGE		
Specificity:	This antibody has been selected as an isotype control as it has no known specificity.		
Source:	Mouse serum		
Molecular Weight:	160 kDa		
Isotype:	IgG		
Extinction Coefficient:	1.36		
Formulation:	Liquid in 0.02 M sodium phosphate, 0.15 M sodium chloride, pH 7.4. Carrier and Preservative Free.		
Purification:	Protein A		
Applications:	ELISA, Western Blot (WB)		
Application Notes:	The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Store at -70 °C for up to one year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.