

## Mouse IgG, Protein A Purified pAb

<b>Catalog No.</b>	CSI20204A CSI20204B	<b>Quantity:</b>	10 mg 50 mg
<b>Description:</b>	Purified from normal serum by immobilized Protein A. >95 percent pure by SDS-PAGE and preservative free.		
<b>Concentration:</b>	10.0 mg/ml		
<b>Specificity:</b>	This antibody has been selected as an isotype control as it has no known specificity.		
<b>Host:</b>	Mouse		
<b>Isotype:</b>	IgG		
<b>Formulation:</b>	Liquid in 0.02M Sodium Phosphate + 0.15M NaCl, pH 7.4		
<b>Purification:</b>	>95 percent pure by SDS-PAGE and preservative free.		
<b>Applications:</b>	ELISA, WB		
<b>Application Notes:</b>	The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	When stored at -80°C, product is stable for 5 years from date of delivery.		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**