cellsciences.com

SERPINE1 Recombinant Human PAI-1 (latent fraction) N-terminal Alexa Fluor 488

Catalog No.	CSI20210A CSI20210B	Quantity:	0.5 mg 1.0 mg
Alternate Names:	PAI, PAI1, PLANH1, SERPINE1		
Description:	Recombinant latent human PAI-1 labeled with Alexa Fluor 488 by iodoacetamide substitution at the N-terminal cysteine		
Concentration:	2.5 mg/ml		
Gene ID:	5054		
Source:	E. coli		
Molecular Weight:	43 kDa		
Formulation:	Frozen liquid in 0.05 M Sodium Phosphate Buffer + 0.1 M NaCl + 1 mM EDTA, pH 6.6.		
Purity:	>95% by SDS-PAGE analysis		
Endotoxin Level:	< 0.1 ng/µg of protein		
Storage & Stability:	Store at -80°C. Protect from light. Stable for 3 years from delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Toll Free: 888-769-1246 Phone: 781-828-0610 Fax: 781-828-0542