

SERPINE1

Recombinant Human PAI-1 (stable mutant) N-terminal Alexa Fluor 594

Catalog No.	CSI20211A CSI20211B	Quantity:	0.5 mg 1.0 mg
Alternate Names:	PAI, PAI1, PLANH1, SERPINE1		
Description:	Recombinant stable mutant human PAI-1 labeled with Alexa Fluor 594 by iodoacetamide substitution at the N-terminal cysteine		
Concentration:	2.5 mg/ml		
Gene ID:	5054		
Source:	<i>E. coli</i>		
Molecular Weight:	43 kDa		
Formulation:	Frozen liquid in 0.05 M Sodium Phosphate Buffer + 0.1 M NaCl + 1 mM EDTA, pH 6.6.		
Purity:	>95% by SDS-PAGE analysis		
Endotoxin Level:	< 0.1 ng/μg of protein		
Storage & Stability:	Store at -80°C. Protect from light. Stable for 3 years from delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.