

## SERPINE1

### Recombinant Human PAI-1 (latent fraction) N-terminal Alexa Fluor 594

<b>Catalog No.</b>	CSI20213A CSI20213B	<b>Quantity:</b>	0.5 mg 1.0 mg
<b>Alternate Names:</b>	PAI, PAI1, PLANH1, SERPINE1		
<b>Description:</b>	Recombinant latent human PAI-1 labeled with Alexa Fluor 594 by iodoacetamide substitution at the N-terminal cysteine		
<b>Concentration:</b>	2.5 mg/ml		
<b>Gene ID:</b>	5054		
<b>Source:</b>	<i>E. coli</i>		
<b>Molecular Weight:</b>	43 kDa		
<b>Formulation:</b>	Frozen liquid in 0.05 M Sodium Phosphate Buffer + 0.1 M NaCl + 1 mM EDTA, pH 6.6.		
<b>Purity:</b>	>95% by SDS-PAGE analysis		
<b>Endotoxin Level:</b>	< 0.1 ng/μg of protein		
<b>Storage &amp; Stability:</b>	Store at -80°C. <b>Protect from light.</b> Stable for 3 years from delivery. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

