

Rabbit IgG Isotype Control Low Endotoxin pAb

Catalog No.	CSI20230A CSI20230B	Quantity:	50 mg 1 g
Description:	Purified from normal serum by immobilized Protein A using low endotoxin methodology and sterile filtered.		
Source:	Rabbit plasma		
Molecular Weight:	160 kDa		
Formulation:	Lyophilized in 0.02M Sodium Phosphate + 0.15M NaCl, pH 7.2		
Purity:	IgG fraction		
Endotoxin Level:	>95% pure by SDS-PAGE and preservative free.		
Target:	This antibody has been selected as an isotype control as it has no known specificity.		
Endotoxin Content:	<10 pg/ul by LAL method		
Reconstitution:	Add deionized water to original volume, aliquot and freeze unused portion.		
Applications:	ELISA WB The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Product is stable for up to 5 years from delivery when stored at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.