

Human Urine Third Trimester (Weeks 29 to 40)

Catalog No.	CSI20367A	Quantity:	50 ml
	CSI20367B		100 ml

Description: Third trimester human Urine collected in week 29-40.
Third trimester human urine collection process:
Human Urine Donor consent forms are obtained and daily monitoring of Human Urine collection is reported. Third trimester human urine is immediately frozen after collection.

Source: Human Urine

Form: Frozen unprocessed liquid.

Storage & Stability: Store at -20°C.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.