

Chicken Anti-Human Alpha-2 Macroglobulin Affinity Purified pAb

| | | | |
|--------------------|-----------|------------------|------|
| Catalog No. | CSI20385A | Quantity: | 1 mg |
| | CSI20385B | | 3 mg |

Description: Polyclonal anti-human IgA procedure: Chickens were immunized with highly purified human serum IgA and eggs were collected. An IgY fraction was prepared and the resulting antibodies were solid phase absorbed to ensure alpha chain specificity. The anti-IgA antibodies were then immunoaffinity purified.

Concentration: Typically greater than 1mg/ml

Purity: Affinity Purified

Host: Chicken

Immunogen: Human IgA

Formulation: Buffer in PBS, pH 7.2.

Purification: Immunoaffinity Purified

Storage & Stability: 2-4°C

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.