

KLK3

Chicken Anti-Human Prostate Specific Antigen Affinity Purified pAb

Catalog No.	CSI20399A CSI20399B	Quantity:	100 µg 0.5 mg
Alternate Names:	Kallikrein-related peptidase 3, Kallikrein-3, Prostate specific antigen, APS, KLK2A1, PSA, hK3		
Description:	<p>Chicken Anti-Human Prostate Specific Antigen (PSA) Polyclonal Antibody.</p> <p>Human Tumor marker Prostate specific antigen (PSA) is a protein produced by the cells of the prostate gland. Prostate Specific Antigen PSA is highly specific for prostatic tissue. Human Prostate Specific antigen PSA is a single chain glycoprotein.</p> <p>Chickens were immunized with purified prostate specific antigen. After multiple immunizations in Freund's adjuvant, eggs were collected and the IgY fraction prepared. The antibodies were then immunoaffinity purified off a prostate specific antigen immunoabsorbant.</p>		
Concentration:	> 1mg/ml		
Gene ID:	354		
Source:	Chicken		
Formulation:	Liquid, PBS pH 7.2		
Purity:	Affinity Purified		
Biological Activity:	The anti-PSA reacts with purified human prostate specific antigen by ELISA and Blotting applications. Optimal working conditions should be determined for your particular assay conditions.		
Storage & Stability:	Store at 2-4°C.		
Certification:	Product manufactured in the USA. All host animals were US donor animals and the human protein was from US origin.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

