

C4, Polyclonal Goat anti-Human pAb, Affinity Purified

Catalog No.	CSI20404A CSI20404B	Quantity:	1 mg 3 mg
Description:	Polyclonal goat anti-Human C4 Procedure: Goat(s) were immunized with highly purified human complement C4 in Freund's adjuvant. The antiserum was solid phase absorbed then the anti-human C4 was immunaffinity purified off an antigen containing immunoabsorbent		
Concentration:	1mg/ml		
Purity:	Affinity Purified		
Specificity:	Polyclonal Antibody Goat anti-Human C4 Specification: This anti-human C4 is monospecific to human C4 by immunoelectrophoresis and gel diffusion techniques		
Host:	Goat		
Immunogen:	Human Compliment C4		
Formulation:	Liquid in PBS, pH 7.2		
Applications:	Western Blot(s) 1:1,000-1:5,000 ELISA 5 ug/ml - 1 ug/ml		
Storage & Stability:	2-8°C		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.