

## Plasminogen Depleted Chicken Serum

<b>Catalog No.</b>	CSI20511A CSI20511B	<b>Quantity:</b>	10 ml 50 ml
<b>Description:</b>	Prepared from frozen chicken serum		
<b>Source:</b>	Chicken Serum		
<b>Purity:</b>	Serum		
<b>Storage &amp; Stability:</b>	When stored at -80°C product is stable for 3 years from date of delivery. <b>Avoid repeated freeze-thaw cycles.</b>		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

