

FGA, FGB

Mouse Fibrinogen Antigen Capture Plate

Catalog No.	CSI20531A CSI20531B	Quantity:	1 plate 5 plates
Alternate Names:	Fib, Fga, Fgb		
Description:	Mouse Fibrinogen Antigen Capture Plate is a 96 well plate with 8 removable strips coated with affinity purified rabbit polyclonal antibody to Mouse Fibrinogen at 100 ul/well and blocked at 300 ul/well. Intended Use: Specific binding of Mouse Fibrinogen antigen for sandwich ELISA experiments.		
Cross-Reactivity:	Cross reacts with Rat Fibrinogen.		
Samples:	Serum, Plasma, Cell culture supernatants		
Sample Size:	100 µl		
Reagents Included:	1 plate. Other reagents not included.		
Storage & Stability:	Store at 2-4°C. Stable for 1 year from delivery.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

