

Plasminogen Depleted Mouse Serum

Catalog No.	CSI20533A CSI20533B	Quantity:	10 ml 50 ml
Description:	Prepared from frozen mouse serum.		
Source:	Mouse Serum		
Purity:	Serum		
Storage & Stability:	When stored at -80°C product is stable for 3 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

