

## Rat Latent PAI-1 Depleted Plasma, Sodium Citrate

<b>Catalog No.</b>	CSI20544A CSI20544B	<b>Quantity:</b>	10 ml 100 ml
<b>Description:</b>	Prepared from frozen rat plasma using immobilized anti-rat PAI-1 IgG.		
<b>Source:</b>	Rat Plasma		
<b>Formulation:</b>	Liquid		
<b>Purity:</b>	Plasma		
<b>Volume:</b>	10.0 ml		
<b>Storage &amp; Stability:</b>	When stored at -80°C product is stable for 3 years from date of delivery. <b>Avoid repeated freeze-thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**