

Rat Renin/Prorenin Double Depleted Plasma

Catalog No.	CSI20545A CSI20545B	Quantity:	10 ml 100 ml
Description:	Prepared from frozen rat plasma using immobilized anti-human Renin IgG.		
Source:	Rat Plasma		
Volume:	1.0 ml		
Formulation:	Liquid in <0.1 ng/ml REN Antigen, <0.1 ng/ml PREN Antigen		
Purity:	Plasma		
Applications:	The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	When stored at -80°C product is stable for 3 years from date of delivery. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

