

Mouse Immunoglobulin G1 Isotype Control Clone MOPC-21 TRI-COLOR

Catalog No. M106 **Quantity:** 0.5 ml

Usage: This isotype control can be used in flow cytometry: Use a maximum of 1 ug per 10^6 cells. Typically 0.5 ug is optimal. Fluorochrome/Protein: 1.

Formulation: 0.5 ml Mouse IgG1 TC in PBS with 50 μ g Ig and EIA-grade BSA and 0.1% sodium azide.
Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.
EIA-grade BSA is added to bring final protein concentration to 4 mg/ml.

Specific Activity: This mouse isotype control can be used for use with single color CD markers.

Storage & Stability: TC conjugates must be stored at 2-8°C...**NEVER STORE FROZEN.**

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

