

Mouse Immunoglobulin G2a Isotype Control Clone 5-205 R-PE

Catalog No. M2A04 **Quantity:** 0.5 ml

Usage: This isotype control can be used in flow cytometry: Use a maximum of 1 ug per 10^6 cells. Typically 0.5 ug is optimal. Fluorochrome/Protein: 1. Excitation 488 nm, peak emission 575 nm.

Formulation: 0.5 ml Mouse IgG2a R-PE in PBS with 50 μ g Ig and EIA-grade BSA and 0.1% sodium azide.
EIA-grade BSA is added to bring final protein concentration to 4-5 mg/ml.
Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.

Specific Activity: This mouse isotype control can be used for use with single color CD markers.

Storage & Stability: R-PE conjugates must be stored at 2-8°C... NEVER STORE FROZEN

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

