

Mouse Anti-Blood Group Antigen H/H Type 2 Clone 19-OLE mAb

Catalog No.	MON2035	Quantity:	100 µg
Concentration:	100 µg/ml.		
Specificity:	The antibody reacts with H type 2 antigen. Source: A BALB/c mouse was immunized with mucin isolated from a mucinous colonic adenocarcinoma. Splenocytes were fused with mouse myeloma SP2/0 cells.		
Host:	Mouse		
Isotype:	IgM		
Clone:	19-OLE		
Formulation:	100 µg purified material in PBS with 0.05% sodium azide. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Applications:	Immunocytochemistry on cultured cells: strongly positive on ethanol fixed epithelial cells. Immunohistochemistry on frozen ethanol fixed section. Immunohistochemistry on paraffin section: very strong positive reaction with ethanol or formalin fixed tissues. MON2035 is also useful in ELISA.		
Storage & Stability:	Store antibody at 2-8°C until expiration date.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

