

## Mouse Anti-Citrate synthetase Clone 8D3-G4 mAb

<b>Catalog No.</b>	MON2057	<b>Quantity:</b>	100 µg
<b>Description:</b>	Mouse Anti-Citrate syntetase Clone 8D3-G4 monoclonal antibody is an enzyme involved in oxidative metabolism (Krebs' cycle) and is abundantly present in the mitochondrial compartment of the cell.		
<b>Alternate Names:</b>	Cis, citrate synthase		
<b>Specificity:</b>	Recognizes citrate synthetase		
<b>Host:</b>	Mouse		
<b>Isotype:</b>	IgG1		
<b>Clone:</b>	8D3-G4		
<b>Formulation:</b>	Lyophilized		
<b>Applications:</b>	Immunohistochemistry (Frozen sections)  The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	Store at 2-4°C until expiration date indicated on the vial. <b>Avoid repeated freeze-thaw cycles.</b>		

**NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.**

