

Anti Chicken IgM + IgG + IgA, clone 408-6.1 Monoclonal Antibody

Catalog No.: MON 5044

Quantity: 1 ml

Specificity

The antibody reacts with anti Chicken IgM + IgG + IgA

Immunoglobulin type

Mouse IgG1

Use

The antibody is suitable for indirect Elisa and immunoblotting.

Instructions for use

The antibody is tested in indirect ELISA with a coating of 50 ng/ml of IgA, IgG and IgM. Also it is tested in competitive ELISA by blocking the binding of mab's to coated Ig's in solution. Can be used in immunoblotting.

Presentation

1 ml with 100 µg purified antibody.

FOR RESEARCH USE ONLY, NOT FOR DRUG, DIAGNOSTIC OR OTHER USE.

