

## Anti Chicken IgM, clone 408-5.1 Monoclonal Antibody

**Catalog No.:** MON 5046

**Quantity:** 1 ml

**Specificity**

The antibody reacts with anti Chicken IgM

**Immunoglobulin type**

Mouse IgG1

**Use**

The antibody is suitable for indirect Elisa and immunoblotting.

**Instructions for use**

The antibody is tested in indirect ELISA with a coating of 50 ng/ml of IgA, IgG and IgM. Also it is tested in competitive ELISA by blocking the binding of mab's to coated Ig's in solution. Can be used in immunoblotting.

**Presentation**

1 ml with 100 µg purified antibody.

FOR RESEARCH USE ONLY, NOT FOR DRUG, DIAGNOSTIC OR OTHER USE.

