

ACRV1

Mouse Anti-Dog Intra-Acrosomal Protein/ACRosomal Vesicle protein 1 Clone Ds-2 mAb

Catalog No.	MON8019	Quantity:	100 µg
Alternate Names:	acrosomal vesicle protein 1		
Description:	This antibody reacts with canine (dog) intra-acrosomal proteins (92 and 42 kDa). Antibody can be used for detection of the state of acrosome of undamaged spermatozoa.		
Concentration:	1 mg/ml		
Gene ID:	607249		
Host:	Mouse		
Isotype:	IgM		
Clone:	Ds-2		
Formulation:	100 µg in solution (PBS containing 15 mM sodium azide, pH 7.4). Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Cross-Reactivity:	Cross reactivity with bovine and pig (pig)		
Applications:	The antibody can be used for immunocytochemistry and Western blotting.		
Application Notes:	For immunocytochemistry membrane permeabilization (acetone) is essential.		
Storage & Stability:	Store at +4°C for six months. For long-term storage aliquot and store at -20°C. Avoid freeze/thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

