

Mouse Anti-Human HL60 mAb Clone IPO-M6 mAb

Catalog No. MON8040 **Quantity:** 100 µg

Description: V17.1 reacts with alkaline phosphatase, an enzyme widely used in immunodetection tests. Source: A BALB/c mouse was immunized with alkaline phosphatase.

Concentration: 100 µg/ml

Host: Mouse

Isotype: IgG1

Clone: IPO-M6

Formulation: 100 µg purified material in PBS with 0.05% sodium azide. **Precaution:** Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.

Applications: Detection of alkaline phosphatase. Amplification of alkaline phosphatase reactions (APAP).

frozen sections	paraffin sections	flow cytometry	ELISA	blocking	Western blotting
+			+		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

