

BD-2, Affinity Purified Human pAb

Catalog No.	PA1362 PA1362BT	Quantity:	100 µg 50 µg
Concentration:	0.5 mg/ml		
Specificity:	Human BD-2		
Host:	Goat		
Immunogen:	Recombinant human BD-2		
Formulation:	Liquid in PBS, pH 7.2, containing 50% glycerol + 1% BSA + 0.02% sodium azide. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Purification:	Affinity chromatography employing immobilized human BD-2 matrix		
Applications:	Western blot: working dilution of 0.5-1.0 µg/ml. The optimal concentration should be determined by the user for each specific application.		
Storage & Stability:	Store at -20°C. For long term storage, aliquot and freeze at -80°C. Avoid repeated freeze-thaw cycles.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

