

CIDEA

Synthetic CIDE-A CT (aa 200-214)

Catalog No. PX016BP **Quantity:** 50 µg

Alternate Names: CIDE-A, cell death activator

Description: Amino acids 200 to 214 of mouse CIDE-A.

This gene encodes the homolog of the mouse protein Cidea that has been shown to activate apoptosis. This activation of apoptosis is inhibited by the DNA fragmentation factor DFF45 but not by caspase inhibitors. Mice that lack functional Cidea have higher metabolic rates, higher lipolysis in brown adipose tissue and higher core body temperatures when subjected to cold. These mice are also resistant to diet-induced obesity and diabetes. This suggests that in mice this gene product plays a role in thermogenesis and lipolysis. Alternatively spliced transcripts have been identified.

Gene ID: 1149

Application: The peptide is used for blocking the activity of anti-mCIDE-A. Incubating the peptide with equal volume of antibody for 30 min at 37°C usually completely blocks the antibody activity in Western blotting.

Formulation: It is supplied as 200 µg/ml, 50 µg/vial, in PBS pH7.2 (10 mM NaH₂PO₄, 10 mM, Na₂HPO₄, 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide.. **Precaution:** Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.

Sequence: GDTEEQSPKPKSTKG

Storage & Stability: Store at -20°C, stable for one year.

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

