

IRAK2

Synthetic Human IL-1 Receptor-Associated Kinase 2 (aa 546-564)(C2) Blocking Peptide

Catalog No.	PX144BP	Quantity:	50 µg
Alternate Names:	IRAK-2, MGC150550, interleukin-1 receptor associated kinase-2		
Description:	IRAK2 encodes the interleukin-1 receptor-associated kinase 2, one of two putative serine/threonine kinases that become associated with the interleukin-1 receptor (IL1R) upon stimulation. IRAK2 is reported to participate in the IL1-induced upregulation of NF-kappaB.		
Gene ID:	3656		
Application:	The peptide is used for blocking the activity of anti-IRAK2. The peptide with equal volume of antibody for 30 min at 37°C usually completely blocks the antibody activity in Western blotting.		
Formulation:	It is supplied as 200 µg/ml, 50 µg/vial , in PBS pH7.2 (10 mM NaH ₂ PO ₄ , 10 mM, Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide.. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Sequence:	APWAGAATPLLPTENGEGR		
Storage & Stability:	Store at -20°C, stable for one year.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

