

TNFRSF17

Synthetic Human BCMA (aa 169-184)(CT) Blocking Peptide

Catalog No.	PX233BP	Quantity:	50 µg
Alternate Names:	BCM, BCMA, CD269, tumor necrosis factor receptor superfamily member 17, B-cell maturation factor, B cell maturation antigen, B-cell maturation protein		
Description:	Amino acids 169 to 183 of human BAFF-R.		
Gene ID:	608		
Application:	The peptide is used for blocking the antibody activity of BCMA. The peptide with equal volume of antibody for 30 min at 37°C usually completely blocks the antibody activity in Western blotting.		
Formulation:	It is supplied as 200 µg/ml, 50 µg/vial , in PBS pH7.2 (10 mM NaH ₂ PO ₄ , 10 mM, Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide. Precaution: Sodium azide is a poisonous and hazardous substance which should be handled by trained staff only.		
Amino Acid Sequence:	PAALSATEIEKSISAR		
Storage & Stability:	Store at -20°C, stable for one year		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

