

Avian Influenza A (H5N1) Hemagglutinin (IN) Peptide

Catalog No.	PX701BP	Quantity: 50 µg
Description:	Location: 14 amino acids near the middle of Avian influenza A (H5N1) hemagglutinin	
Formulation:	It is supplied as 50 µg at 200 µg/ml in PBS pH7.2 (10 mM NaH ₂ PO ₄ , 10 mM Na ₂ HPO ₄ , 130 mM NaCl) containing 0.1% bovine serum albumin and 0.02% sodium azide.	
Storage & Stability:	Store at -20°C, stable for one year.	
Applications:	The peptide is used for blocking the activity of the anti-avian influenza A (H5N1) hemagglutinin antibody (catalog number PX701A). It usually blocks the antibody activity completely in ELISA assays by incubating the peptide with equal volume of antibody for 30 min at 37°C.	

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.



Cell Sciences, Inc.
480 Neponset Street
Bldg 12A
Canton, MA 02021

Toll Free: 888-769-1246
Phone: 781-828-0610
Fax: 781-828-0542

E-mail: info@cellsciences.com
Web site: www.cellsciences.com