

## Chicken Anti-human GDNF pAb

<b>Catalog No.</b>	PX739	<b>Quantity:</b>	0.5 mg
<b>Alternate Names:</b>	Glial cell line-derived neurotrophic factor, Astrocyte-derived trophic factor, hGDNF, ATF1, ATF2, HSCR3, HFB1-GDNF		
<b>Description:</b>	Chicken Anti-human GDNF polyclonal antibody		
<b>Specificity:</b>	Human GDNF		
<b>Host:</b>	Chicken		
<b>Immunogen:</b>	GDNF mixed peptides		
<b>Cross-Reactivity:</b>	Mouse, Rat		
<b>Applications:</b>	ELISA, Western blot		
<b>Application Notes:</b>	The optimal concentration should be determined by the user for each specific application.		
<b>Storage &amp; Stability:</b>	Store at 2-8 °C for up to one week or at -20 °C as an undiluted liquid for up to two years.		

NOT FOR HUMAN USE. FOR RESEARCH ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.